

## Preparation of de yolked embryos

Ver.2019-3-26

Zebrafish embryos (500 – 2000) in 3.5 cm dish

- ↓ reduce the culture medium to about 1 mL
- ↓ add 1 mL Pronase stock
- ↓ 28.5°C, 15 min, in incubator
- ↓ wash with ice cold E3 medium, thoroughly
- ↓ transfer to another dish with blue tip (for de-chorion)
- ↓ wash with ice cold deyolking buffer, 3 times
- ↓ deyolk by pipetting with yellow tip (check with microscope)
- ↓ transfer to Eppen tube(s) (two tubes for >1000 embryos)
- ↓ fill up with ice cold deyolking buffer
- ↓ mix by shaker (EYELA, CM-1000), 5 min, speed 8
- ↓ cfg., 400 x g, 4°C, 30 sec

ppt.

- ↓ suspend in 1.5 mL ice cold wash buffer
- ↓ mix by shaker, 5 min, speed 8
- ↓ cfg., 400 x g, 4°C, 30 sec

ppt.

- ↓ suspend in 1.5 mL ice cold PBS
- ↓ mix by shaker, 5 min, speed 8
- ↓ cfg., 400 x g, 4°C, 30 sec

ppt.

- ↓ store at -70°C

### ➤ Pronase stock (for zebrafish embryo)

#### • Stock solution

5 mg/mL Pronase, 0.05% Phenol red in E3 medium

#### • Preparation

10 mg Phenol red

↓ 20 mL E3 medium

↓ add several drops of 1 mM NaHCO<sub>3</sub> to turn red color

↓ dissolve 100 mg Pronase

↓ aliquot each 1 mL in Eppen tubes

↓ store at -20°C